

# PGT Digital Imaging...

**...ImageCraft™, highest pixel density per displayed image... over 20% more than its competition.**

At the touch of a key... micrographs digitized, color enhanced, and processed to provide more data than you ever thought possible.

PGT **ImageCraft™**... 400 x 400 pixel image display with immediate access to 512 x 400 pixel images.

**MegaPlane™**... ultra high resolution imaging display system... features immediate access to 1024 x 800 pixel images, 256 simultaneously assigned colors from 16 million color palette... including 256 true gray levels on monitor and *instantaneous* video refresh.

Image contrast enhancement, automatic area fraction display, instantaneous zoom and roam, sophisticated image comparison capabilities... and many special purpose software options available for a variety of image analysis requirements.

There is a PGT system especially suited to your needs... with EDS X-ray analysis functions, or specifically for image analysis.

Remember these names...  
**SYSTEM4plus - ImageCraft™**  
**- Ultra - MegaPlane™**



**PRINCETON GAMMA-TECH, INC.**

1200 State Road, Princeton, NJ 08540 • 609-924-7310 • Telex: (WUT) 843486 GAMMATECH PRIN, (ITT) 4754029 PGT USA  
568 Weddell Drive, Suite 1, Sunnyvale, CA 94089 • 408-734-8124

# CamScan

## Scanning Electron Microscopes

Technical Excellence, Outstanding Quality and Reliable Service.



The CamScan line of advanced research scanning microscopes is the result of extensive programs of design, development and engineering based on a real appreciation of the needs of the microscopist. Nothing has been spared to obtain the highest technical specifications and workmanship.

### CamScan Key Features:

#### LARGE SPECIMEN CHAMBER


Large enough to house virtually any sample encountered in scanning microscopy. A viewing window and internal illuminator simplify the manipulation of bulky or awkwardly shaped samples.

#### EUCENTRIC STAGE

The technical superiority of the eucentric stage design provides high stability precision movements with excellent control of position and orientation for specimens of all types and sizes.

#### MICROANALYTICAL GEOMETRY

Designed specifically to give the optimum geometrical configuration between sample, electron beam and X-ray detectors for accurate quantitative microanalysis.



#### UNIQUE ERGONOMIC DESIGN

The layout and functions of the microscope controls ensure optimum performance, rapid results and simple operation for all samples.



**WEST**  
6678 Owens Drive  
Pleasanton, California 94566  
(415) 463-0200 • Tlx 704483

**EAST**  
780 Commonwealth Drive  
Warrendale, Pennsylvania 15086  
(412) 776-5260 • Tlx 902867

For more information on CamScan Scanning Electron Microscopes, please write for our brochure.

# PUBLICATIONS FROM 1985 MRS FALL MEETING

The following prices are effective until **March 30, 1986**. Prices after this date will be higher.  
Order your copies today to receive priority shipment.

PROCEEDINGS		MRS Members	U.S.	Foreign	Quantity	Total
		& Fall Meeting Attendees	Nonmembers	Nonmembers		
A	Beam-Solid Interactions and Phase Transformations	\$32	\$38	\$43	_____	_____
B	Rapid Thermal Processing	\$30	\$36	\$41	_____	_____
C	Semiconductor on Insulator and Thin Film Transistor Technology	\$32	\$38	\$43	_____	_____
E	Thin Films—Interfaces and Phenomena	\$39	\$47	\$52	_____	_____
G	Biomedical Materials	\$33	\$40	\$45	_____	_____
H	Layered Structures and Epitaxy	\$34	\$41	\$46	_____	_____
I	Phase Transitions in Condensed Systems—Experiments and Theory	\$32	\$38	\$43	_____	_____
J	Rapidly Solidified Alloys and Their Mechanical and Magnetic Properties	\$28	\$34	\$39	_____	_____
K	Oxygen, Carbon, Hydrogen, and Nitrogen in Crystalline Silicon	\$32	\$38	\$43	_____	_____
L	Defect Properties and Processing of High-Technology Nonmetallic Materials	\$33	\$40	\$45	_____	_____
P	Defects in Glasses	\$28	\$34	\$39	_____	_____
Q	Materials Problem Solving with the Transmission Electron Microscope	\$33	\$40	\$45	_____	_____
R	Computer-Based Microscopic Description of the Structure and Properties of Materials	\$30	\$36	\$41	_____	_____
S	Cement-Based Composites: Strain Rate Effects on Fracture	\$24	\$30	\$35	_____	_____
T	Fly Ash and Coal Conversion By-Products: Characterization, Utilization and Disposal II	\$24	\$30	\$35	_____	_____
Y	Frontiers in Materials Education	\$24	\$30	\$35	_____	_____
<b>EXTENDED ABSTRACTS</b>						
D	Beam Induced Chemical Processes	\$5	\$8	\$8	_____	_____
N	Fractal Aspects of Materials	\$5	\$8	\$8	_____	_____
O	Nonlinear Optical Materials	\$5	\$8	\$8	_____	_____
					<b>TOTAL</b>	_____

*Prepayment required.*

Return to: Publications Department, Materials Research Society, 9800 McKnight Road, Suite 327, Pittsburgh, Pa 15237;  
telephone (412) 367-3012



**MATERIALS RESEARCH SOCIETY**

9800 McKnight Road, Suite 327  
Pittsburgh, Pennsylvania 15237  
U.S.A.

JAN 2 1985

NON-PROFIT ORG.  
U.S. Postage  
**PAID**  
Pittsburgh, PA  
Permit No. 1807

JOHN B BALLANCE ( )  
MATERIALS RESEARCH SOCIETY  
9800 MCKNIGHT ROAD  
SUITE 327  
PITTSBURGH, PA 15237  
USA

## COMING IN THE NEXT ISSUE

- MRS in 1985—The Year in Retrospect
- Summary of 1985 MRS Fall Meeting Symposia
- Preview of 1986 MRS Spring Meeting
- Details on 1986 Short Courses  
(to be held in conjunction with 1986 Spring Meeting)